Regensburg, den 5. Januar 2004

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Dear Participant, dear Colleague,

Please find enclosed a brief analysis report of the recent proficiency test panel "Bacterial Genome Detection (PCR / NAT)". It summarizes some benchmarks and the overall assessment of results reported by all of the the participating laboratories.

The accuracy of your individual results is indicated on your personal certificate issued by INSTAND e.V. (attached to this letter).

This highly desired program for external validation was activated by the *German Society of Hygiene and Microbiology* (DGHM) and is now organized by INSTAND e.V., Düsseldorf, Germany. As you can learn from "www.instand-ev.de", this new segment of their huge proficiency testing program is now open for diagnostic laboratories worldwide. We are aiming at two validation rounds per year (April / November), a reasonable price, and a permanently expanding coverage of relevant bacterial pathogens. Any general or specific comments on sample composition or suggestions on additional bacterial, fungal, or parasitic organisms which should be covered, are highly appreciated.

Please do not feel irritated by some strong-positive samples among the present set. To make participants and ourselves familiar with the concept and the suitability of the proprietary sample matrix, "passing" should be easy in the initial rounds of proficiency testing. But as the program continues, it is our challenge (and duty) to rise at least some of the hurdles with every round of external validation.

Further information on this proficiency test panel "Bacterial Genome Detection (PCR / NAT)" can be downloaded in form of pdf-files at "www.udo-reischl.de", subsection "INSTAND Ringversuche (PCR / NAT)". Although the preferred language of these documents is German, we are aiming to provide at least some key documents and the tables in a bilingual style.

Hopefully you found this kind of external validation benefitting and you will continue to join the future rounds of our QC program in bacterial genome detection (NAT/ PCR).

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Sample composition and expected results.

Gruppe A	Erwartet / expected		Probenzusammensetzung / Sample composition
32301	∅	62	<i>Escherichia coli</i> K12
32302	+	61 / 72	<i>Helicobacter pylori</i> (~ 1x10 ⁴ CFU/mL); Clarithromycin susceptible (wildtype 28S rDNA sequence)
32303	∅	62	<i>Helicobacter mustelae</i> (~ 1x10 ⁵ CFU/mL);
32304	++	61 / 71	<i>Helicobacter pylori</i> (~ 1x10 ⁵ CFU/mL) Clarithromycin resistant (GAG mutation in 28S rDNA)

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Absolute numbers of reported individual results.

n = 11	Probennummer (Sample no.)					Inhibition			
	32301	32302	32303	32304		32301	32302	32303	32304
Befund Result									
Positiv	0	6	7	10	n.d.	0	0	0	0
Negativ	10	4	3	0	nein no	11	11	11	11
Fraglich Questionable	1	1	1	1	ja yes	0	0	0	0

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Absolute numbers and relative frequency of reported true positive and true negative results among various NAT methods.

NAT-Methode [Code] (total number *)	NAT richtig positiv True positive results			NAT richtig negativ True negative results		
	Absolut Absolute	Relativ Relative	%	Absolut Absolute	Relativ Relative	%
In house PCR assay [28] (n = 11)	18	18 / 22	82	15	15 / 22	68 ¹⁾
Andere / other [29] (n = 0)	-	-	-	-	-	-

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Comments:

¹⁾ Seven false-positive results were observed for sample #31303 (*H. mustelae*);
 These participants have either used a ribosomal gene or sections within the urease gene as target sequence for their "*H. pylori*-specific" PCR assays.

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Results reported by our Reference Laboratories

Probenversand / shipment of samples: 12. November 2003

(results of the QC-program reference labs; listed but not depicted in alphabetical order):

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status 11.2003

Reischl / Lehn / Wolf



Evaluation (quantitative Real-Time PCR):

- 1 Hepy 32301/neu
- 2 Hepy 32302/neu >41.00
- 3 Hepy 32303/neu
- 4 Hepy 32304/neu 28.46
- 9 poKo Hepy/Urease 18.73
- 10 neg. MM H2O

LightCycler PCR protocol:

Reischl, U., B. Leppmeier, M. Heep, D. Beck, and N. Lehn (2001) Rapid and specific detection of *Helicobacter pylori* by LightCycler PCR. In: Rapid Cycle Real-Time PCR: Methods and Applications (Meurer, S., Wittwer, C., and Nakagawra, K.-I., eds.), ISBN 3-540-66736-9, Springer Press, Heidelberg, pp. 323-330.

U. Reischl/RIMMH/11.2003

433 Bakteriengenom-Nachweis *Helicobacter pylori*

status 11.2003

Reischl / Lehn / Wolf



Evaluation (quantitative Real-Time PCR):

- 7 Hepy 32301
- 8 Hepy 32302 33.60
- 9 Hepy 32303 25.98
- 10 Hepy 32304 23.89
- 11 Pos Hepy 16S 15.28
- 12 Neg Ko


LightCycler PCR protocol:

Reischl, U., B. Leppmeier, M. Heep, D. Beck, and N. Lehn (2001) Rapid and specific detection of *Helicobacter pylori* by LightCycler PCR. In: Rapid Cycle Real-Time PCR: Methods and Applications (Meurer, S., Wittwer, C., and Nakagawra, K.-I., eds.), ISBN 3-540-66736-9, Springer Press, Heidelberg, pp. 323-330.

U. Reischl/RIMMH/11.2003

433 Bakteriengenom-Nachweis *Helicobacter pylori*

status 11.2003

Evaluation (quantitative Real-Time PCR):


Reischl / Lehn / Wolf

H. pylori
 Clari-Resistance
 28S rDNA

- 1 Hepy 32301
- 2 Hepy 32302
- 3 Hepy 32303
- 4 Hepy 32304
- 5 Plasmid AA
- 6 Plasmid AG
- 7 Plasmid GA
- 8 Plasmid CA
- 9 neg. MM H2O

LightCycler PCR protocol:
 unpublished *in house* protocol.

- Centrifuge the vial to pellet the lyophilized sample and add 300 µL of PCR-grade water.
- Incubate at room temperature overnight or at 37°C for 30 min.
- Vortex and process an aliquot of each sample according to your sample preparation and diagnostic PCR / NAA protocol.

U. Reischl/RIMMH/11.2003


432- 438 Bakteriengenom-Nachweis *all sets*

status 11.2003

➤ Evaluation (quantitative Real-Time PCR):

Reischl / Lehn / Wolf

1	RV-Bope 32201	27.68
2	RV-Bope 32202	27.04
3	RV-Bope 32203	26.95
4	RV-Bope 32204	27.71
5	RV-Hepy 32301	28.17
6	RV-Hepy 32302	29.03
7	RV-Hepy 32303	28.84
8	RV-Hepy 32304	25.88
9	RV-EHEC 32401	28.74
10	RV-EHEC 32402	28.85
11	RV-EHEC 32403	28.70
12	RV-EHEC 32404	27.72
13	RV-Bobu 32501	27.84
14	RV-Bobu 32502	28.04
15	RV-Bobu 32503	28.14
16	RV-Bobu 32504	28.43
17	RV-Lepn 32601	27.85
18	RV-Lepn 32602	27.58
19	RV-Lepn 32603	28.29
20	RV-Lepn 32604	27.79
21	RV-Salspp. 32701	28.98
22	RV-Salspp. 32702	28.82
23	RV-Salspp. 32703	27.88
24	RV-Salspp. 32704	28.08
25	RV-Lisp 32801	27.97
26	RV-Lisp 32802	24.56
27	RV-Lisp 32803	28.10
28	RV-Lisp 32804	28.50
29	poKo β-Globin 75ng	23.12
30	neg. MM H2O	37.42

LightCycler PCR protocol:
 LightCycler Control Kit DNA
 Roche Cat. No. 2 158 833

human cells / PCR reaction: ~10³

030923_3_RV_β-Globin-ampl.bmp

- Centrifuge the vial to pellet the lyophilized sample and add 300 µL of PCR-grade water.
- Incubate at room temperature overnight or at 37°C for 30 min.
- Vortex and process an aliquot of each sample according to your sample preparation and diagnostic PCR / NAA protocol.

U. Reischl/RIMMH/11.2003


432- 438 Bakteriengenom-Nachweis *all sets*

status 11.2003

➤ Evaluation (quantitative Real-Time PCR):

Reischl / Lehn / Wolf

1	RV-Bope 32201	23.73
2	RV-Bope 32202	23.57
3	RV-Bope 32203	24.50
4	RV-Bope 32204	23.50
5	RV-Hepy 32301	23.87
6	RV-Hepy 32302	24.30
7	RV-Hepy 32303	24.21
8	RV-Hepy 32304	24.10
9	RV-EHEC 32401	23.05
10	RV-EHEC 32402	24.88
11	RV-EHEC 32403	22.33
12	RV-EHEC 32404	23.07
13	RV-Bobu 32501	23.66
14	RV-Bobu 32502	23.50
15	RV-Bobu 32503	23.81
16	RV-Bobu 32504	23.90
17	RV Lepn 32601	23.56
18	RV Lepn 32602	23.45
19	RV Lepn 32603	24.24
20	RV Lepn 32604	23.65
21	RV Salspp. 32701	23.32
22	RV Salspp. 32702	23.60
23	RV Salspp. 32703	23.62
24	RV Salspp. 32704	24.47
25	RV Lisp 32801	23.75
26	RV Lisp 32802	22.72
27	RV Lisp 32803	23.68
28	RV Lisp 32804	24.32
29	poKo parC (E.coli)	18.76
30	neg. MM H2O	


LightCycler PCR protocol:
 -unpublished-

Enterobact. / PCR reaction: ~10⁵

030924_3_RV_parC-ampl.bmp

- Centrifuge the vial to pellet the lyophilized sample and add 300 µL of PCR-grade water.
- Incubate at room temperature overnight or at 37°C for 30 min.
- Vortex and process an aliquot of each sample according to your sample preparation and diagnostic PCR / NAA protocol.

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Detection of *Helicobacter pylori* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)

Instructions for testing and Code numbers (2 pages)

Precautions

The samples in this panel were subjected to inactivation and they are **not** considered hazardous or capable of transmitting infection.

Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *H. pylori*-PCR/NAA established in your routine diagnostic setting.

Testing and reporting of results

The concept of this proficiency panel is designed for determination of *H. pylori* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like β-globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers.

Although only reporting of qualitative results is requested, you are free to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

Completing of Report Form

Please return the completed Report Form no later than the requested date (postmark date) to:

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BACTERIAL GENOME DETECTION

CODE Numbers PCR-/NAA *Helicobacter pylori* (433)

Code numbers are subdivided into seven groups. Although providing technical specifications of the applied assays is highly desired for statistical purposes (individual data will be kept confidential), completing groups [I] to [V], and [VII], is voluntary.

Completing of group [VI] is mandatory.

Group [I] (Sample preparation)

11: Proteinase K / spin column
12: Commercial DNA extraction kit **
13: Phenol / chloroform extraction 19: Other **

Group [II] (Amplification)

27: Commercial assay system / kit **
28: *In house* PCR assay 29: Other **

Group [III] (Detection / identification)

31: Commercial assay system (codes 21-27)
32: Agarose gel electrophoresis
33: Hybridization with labelled probe
34: Nested PCR
35: Real-Time PCR (TaqMan format)
36: Real-Time PCR (LightCycler)
37: DNA sequencing 39: Other **

Group [IV] (Positive and/or inhibition control)

41: Commercial assay system (codes 21-27)
42: Internal control (recombinant plasmid etc.)
43: External controls (plasmid, β -globin gene, etc.)
44: Genomic DNA of the corresponding pathogen
48: Other ** 49: No inhibition control applied

Group [V] (Target gene)

51: Commercial assay system (codes 21-27)
52: Ribosomal gene (16S rDNA, 28S rDNA, ITS)
53: Urease gene 59: Other **

Group [VI] (Results)

61: **Positive** 63: **Questionable**
62: **Negative** 64: **Inhibition**

Gruppe [VII] (Molecular susceptibility testing)

71: Presumably Clarithromycin-resistant
72: Presumably Clarithromycin-susceptible

** (please specify in the **Comments** section)