# **INSTAND**

## **BEGLEITHEFT**

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## BACTERIAL GENOME DETECTION

### PCR-/NAT *Chlamydia trachomatis* & *Neisseria gonorrhoeae* (RV 530)

#### Detection of *Chlamydia trachomatis* & *Neisseria gonorrhoeae* by PCR or other procedures for nucleic acid amplification & detection (NAA)

#### Instructions for testing and Code numbers (2 pages)

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *C. trachomatis* & *N. gonorrhoeae*-PCR/NAA established in your routine diagnostic setting.

Note: when using commercial assay systems, 100 µl aliquots of the renaturated samples can be processed like the sediment of centrifuged urine specimens or filled up with sterile water to reach the minimal input volume requested by the workup scheme of some commercial assays.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *C. trachomatis* and *N. gonorrhoeae* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

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Although only reporting of qualitative results is requested, you are free to provide quantitative or semiquantitative data (e.g., +, ++, etc.) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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#### Group [II] (Amplification)

20: GenProbe CT/NG    21: LightMix CT/NG (TIB Molbiol)  
22: Roche COBAS      23: Xpert CT/NG (Cepheid)  
24: BD ProbeTec      25: artus CT  
26: Abbott RealTime CT/NG  
27: Other commercial assay / kit \*\*  
28: *In house* PCR assay      29: Other \*\*

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31: Commercial assay system (codes 20-27)  
32: Agarose gel electrophoresis  
33: Hybridization with labelled probe  
34: Nested PCR  
35: Real-Time PCR  
37: DNA sequencing      39: Other \*\*

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#### Group [VI] (Results\*)

61: Positive Ct      64: Negative  
62: Positive Ct and GO      65: Questionable  
63: Positive GO      66: Inhibition

\* Ct: *C. trachomatis*; GO: *N. gonorrhoeae*  
\*\* (please specify in the **Comments** section)

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
## BACTERIAL GENOME DETECTION

### PCR-/NAA *Chlamydia trachomatis* (RV 531)

#### Detection of *Chlamydia trachomatis* by PCR or other procedures for nucleic acid amplification & detection (NAA)

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *C. trachomatis*-PCR/NAA established in your routine diagnostic setting.

Note: when using commercial assay systems, 100 µl aliquots of the renatured samples can be processed like the sediment of centrifuged urine specimens.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *C. trachomatis* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

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
## BACTERIAL GENOME DETECTION

### PCR-/NAA *Bordetella pertussis* (RV 532)

**Detection of *Bordetella pertussis* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *B. pertussis*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *B. pertussis* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers.

Although only reporting of qualitative results is requested, you are free to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

#### Completing of Report Form

Please return the completed Results Form no later than the requested date (postmark date) to:

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#### Group [V] (Target gene)

51: Commercial assay system (codes 20-27)  
52: Ribosomal gene (16S rDNA, 28S rDNA, ITS)  
53: Insertion sequence *IS* 481  
54: Pertussis toxin gene      59: Other \*\*

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## BACTERIAL GENOME DETECTION

### PCR-/NAA *Helicobacter pylori* (RV 533)

**Detection of *Helicobacter pylori* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *H. pylori*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *H. pylori* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

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
## BACTERIAL GENOME DETECTION

### PCR-/NAA EHEC / STEC (RV 534)

#### Detection of EHEC / STEC DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and EHEC/STEC-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of EHEC/STEC DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

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## BACTERIAL GENOME DETECTION

### PCR-/NAA *Borrelia burgdorferi* (RV 535)

**Detection of *Borrelia burgdorferi* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *B. burgdorferi*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *B. burgdorferi* sensu lato DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

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## BACTERIAL GENOME DETECTION

### PCR-/NAA *Legionella pneumophila* (RV 536)

**Detection of *Legionella pneumophila* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *L. pneumophila*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *L. pneumophila* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like β-globin, etc.).

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## BACTERIAL GENOME DETECTION

### PCR-/NAA *Salmonella enterica* (RV 537)

**Detection of *Salmonella enterica* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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When the pellets are completely dissolved, the resulting suspensions should be considered as enrichment broths or native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *Salmonella enterica*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *Salmonella enterica* DNA in the sample material. RNA should be detectable as well. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like β-globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of *S. enterica* is requested, you are free to report the corresponding serovars and/or to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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#### Group [II] (Amplification)

20: RealArt *Salmonella* kit (Artus)  
21: foodproof *Salmonella* kit (Biotecon)  
27: Commercial assay system / kit \*\*  
28: *In house* PCR assay      29: Other \*\*

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#### Group [V] (Target gene)

51: Commercial assay system (codes 20-27)  
52: Ribosomal gene (16S rDNA, 28S rDNA)  
53: *himA* gene      54: *gyrB* gene      55: *invA* gene  
56: Ribosomal ITS region      59: Other \*\*

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
## BACTERIAL GENOME DETECTION

### PCR-/NAA *Listeria spp.* (RV 538)

**Detection of *Listeria spp.* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *Listeria spp.* - PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *Listeria spp.* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like β-globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of enrichment broths or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of *Listeria spp.* is requested, you are free to report the corresponding species names and/or to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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
## BACTERIAL GENOME DETECTION

### PCR-/NAA MRSA / CA-MRSA (RV 539)

**Direct detection of MRSA / CA-MRSA DNA**  
by PCR or other procedures for nucleic  
acid amplification & detection (NAA)

**Instructions for testing and**  
Code numbers (2 pages)

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#### Specimen preparation

Samples should be stored at 4°C upon receipt.

The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing.

When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and MRSA / CA-MRSA PCR/NAA assays established in your routine diagnostic setting. **Note:** PVL-(Panton-Valentine leukocidin) positive MRSA isolates are termed **CA-**(community acquired-) **MRSA**.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of MRSA and/or CA-MRSA DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of enrichment broths or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of MRSA is requested, you are free to report the corresponding species names and/or the presence of the *lukFS* gene (CA-MRSA) [Code 71] in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

#### Completing of Report Form

Please return the completed Report Form no later than the requested date (postmark date) to:  
INSTAND e.V.

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#### Group [II] (Amplification)

20: BD MAX / BD GeneOhm MRSA (Becton Dickinson)  
21: FT MRSA (Hain Lifescience)  
22: RIDAGENE MRSA (r-Biopharm)      23: -  
24: Xpert MRSA / GeneXpert (Cepheid)  
25: LightCycler MRSA Advanced (Roche)  
26: LightMix MRSA (TIB Molbiol)  
27: Commercial assay system / kit \*\*  
28: *In house* PCR assay      29: Other \*\*

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31: Commercial assay system (codes 20-27)  
32: Agarose gel electrophoresis  
33: Hybridization with labelled probe  
34: Nested PCR,      36: Real-Time PCR (LightCycler)  
37: DNA sequencing      39: Other \*\*

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## BACTERIAL GENOME DETECTION

### PCR-/NAA *Chlamydia pneumoniae* (RV 540)

**Detection of *C. pneumoniae* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **1000 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and corresponding aliquots processed as such using typical protocols for DNA extraction and *Chlamydia pneumoniae*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *Chlamydia pneumoniae* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like β-globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of enrichment broths or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of *C. pneumoniae* is requested, you are free to report the corresponding species names and/or to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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## BACTERIAL GENOME DETECTION

### PCR-/NAA *Mycoplasma pneumoniae* (RV 541)

**Detection of *M. pneumoniae* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and corresponding aliquots processed as such using typical protocols for DNA extraction and *M. pneumoniae*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *Mycoplasma pneumoniae* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like β-globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of enrichment broths or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of *M. pneumoniae* is requested, you are free to report the corresponding species names and/or to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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## BACTERIAL GENOME DETECTION

### PCR-/NAA *Coxiella burnetii* & *Bacillus anthracis* (RV 542)

**Detection of *C. burnetii* and *B. anthracis* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and corresponding aliquots processed as such using typical protocols for DNA extraction and *C. burnetii*- and/or *B. anthracis*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *Coxiella burnetii* & *B. anthracis* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of enrichment broths or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of *Coxiella burnetii* and *B. anthracis* is requested, you are free to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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- 11: Proteinase K / spin column
- 12: Commercial DNA extraction kit \*\*
- 13: Phenol / chloroform extraction      19: Other \*\*

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**Group [III] (Detection / identification)**

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- 31: Commercial assay system (codes 20-27)
- 32: Agarose gel electrophoresis
- 33: Hybridization with labelled probe
- 34: Nested PCR
- 35: Real-Time PCR
- 37: DNA sequencing      39: Other \*\*

**Gruppe [IV] (Positiv u/o Inhibitionskontrollen)**

**Group [IV] (Positive and/or inhibition control)**

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- 41: Commercial assay system (codes 20-27)
- 42: Internal control (recombinant plasmid etc.)
- 43: External controls (plasmid,  $\beta$ -globin gene, etc.)
- 44: Genomic DNA of the corresponding pathogen
- 48: Other \*\*      49: No inhibition control applied

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**Group [V] (Target gene)**

- 51: Kommerzielles Testsystem (Codes 20-27)
- 52: Ribosomales Gen (16S rDNA, 23S rDNA)
- 53: Transposase Gen (IS1111)    54: COM-1 Gen
- 55: rpoB Gen    56: pX01    57: pX02    59: Andere \*\*

- 51: Commercial assay system (codes 20-27)
- 52: Ribosomal gene (16S rDNA, 23S rDNA)
- 53: Transposase gene (IS1111)    54: COM-1 gene
- 55: rpoB gene    56: pX01    57: pX02    59: Other \*\*

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**Group [VI] (Results)**

- 61: **Positiv *C. burnetii***      64: **Negativ**
- 62: **Positiv Beide**      65: **Fraglich**
- 63: **Positiv *B. anthracis***    66: **Inhibition**

- 61: **Positive *C. burnetii***      64: **Negative**
- 62: **Positive Both**      65: **Questionable**
- 63: **Positive *B. anthracis***    66: **Inhibition**

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optional: 71: only *C. burnetii* assay performed  
72: only *B. anthracis* assay performed

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

## BACTERIAL GENOME DETECTION

### PCR-/NAA *Francisella tularensis* (RV 543)

**Detection of *F. tularensis* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)**

**Instructions for testing and Code numbers (2 pages)**

#### Precautions

The samples in this panel were subjected to inactivation and they are **not** considered hazardous or capable of transmitting infection.

#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and corresponding aliquots processed as such using typical protocols for DNA extraction and *F. tularensis*-PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *Francisella tularensis* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like β-globin, etc.).

Please apply the same methods for DNA extraction and amplification like you are using for the workup of enrichment broths or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of *F. tularensis* is requested, you are free to report the corresponding species names and/or to provide quantitative or semiquantitative data (e.g., +, +++) in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

#### Completing of Report Form

Please return the completed Report Form no later than the requested date (postmark date) to:

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## BACTERIAL GENOME DETECTION

### CODE Numbers PCR-/NAA *F. tularensis* (543)

Code numbers are subdivided into six groups. Although providing technical specifications of the applied assays is highly desired for statistical purposes (individual data will be kept confidential), completing groups [I] to [V] is on a voluntary basis.

Completing of group [VI] is mandatory.

#### Group [I] (Sample preparation)

11: Proteinase K / spin column  
12: Commercial DNA extraction kit \*\*  
13: Phenol / chloroform extraction      19: Other \*\*

#### Group [II] (Amplification)

20: LightMix *Francisella tularensis* kit (TIB Molbiol)  
  
27: Commercial assay system / kit \*\*  
28: *In house* PCR assay      29: Other \*\*

#### Group [III] (Detection / identification)

31: Commercial assay system (codes 20-27)  
32: Agarose gel electrophoresis  
33: Hybridization with labelled probe  
34: Nested PCR  
35: Real-Time PCR (TaqMan format)  
36: Real-Time PCR (LightCycler)  
37: DNA sequencing      39: Other \*\*

#### Group [IV] (Positive and/or inhibition control)

41: Commercial assay system (codes 20-27)  
42: Internal control (recombinant plasmid etc.)  
43: External controls (plasmid,  $\beta$ -globin gene, etc.)  
44: Genomic DNA of the corresponding pathogen  
48: Other \*\*      49: No inhibition control applied

#### Group [V] (Target gene)

51: Commercial assay system (codes 20-27)  
52: Ribosomal gene (16S rDNA, 23S rDNA)  
53: *tul4/lpnA* gene  
55: *fopA* gene      59: Other \*\*

#### Group [VI] (Results)

61: Positive      63: Questionable  
62: Negative      64: Inhibition

\*\* (please specify in the **Comments** section)

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
## BACTERIAL GENOME DETECTION

### PCR-/NAA Carbapenemase genes (RV 544)

Detection of Carbapenemase-encoding genes in enterobacteriaceae by PCR or other procedures for nucleic acid amplification & detection (NAA)

Instructions for testing and Code numbers (2 pages)

#### Precautions

The samples in this panel were subjected to inactivation and they are **not** considered hazardous or capable of transmitting infection.

#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as culture aliquots or native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and carbapenemase gene PCR/NAA assays established in your routine diagnostic setting. **Note:** The DNA content of the samples may be significantly lower compared to direct extraction or boiling a bacterial colony.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of DNA from bacterial carbapenemase genes in the samples. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to target organisms, thereby allowing inhibition controls targeting human genes like  $\beta$ -globin. Please apply the same methods for DNA extraction and amplification like you are using for the workup of cultures or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of carbapenemase genes is requested, you are free to report the corresponding carbapenemase genes [codes 70 to 79] in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

#### Completing of Report Form

Please return the completed Report Form no later than the requested date (postmark date) to:  
INSTAND e.V.  
Abt. Qualitätssicherung  
Postfach 250211  
40093 Düsseldorf, GERMANY.

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#### Gruppe [II] (Amplifikation)

20: Check MDR Carba test (Fa. Checkpoints, distributed by HAIN Lifescience)  
21: Check Direct CPE *oder* MDR (Fa. Checkpoints, distributed by HAIN Lifescience)  
22: hyplex Superbug      23: hyplex Eazyplex  
24: LightMix (TIB Molbiol)  
27: Kommerzielles Testsystem / Kit \*\*  
28: *In house* PCR Protokoll      29: Andere \*\*

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## BACTERIAL GENOME DETECTION

### CODE Numbers PCR-/NAA Carbapenemases (544)

Code numbers are subdivided into six groups. Although providing technical specifications of the applied assays is highly desired for statistical purposes (individual data will be kept confidential), completing groups [I] to [V] is on a voluntary basis.

Completing of group [VI] is mandatory.

#### Group [I] (Sample preparation)

11: Proteinase K / spin column  
12: Commercial DNA extraction kit \*\*  
13: Phenol / chloroform extraction      19: Other \*\*

#### Group [II] (Amplification)

20: Check MDR Carba test (Fa. Checkpoints)  
21: Check Direct CPE *or* MDR (Fa. Checkpoints)  
22: hyplex Superbug      23: hyplex Eazyplex  
24: LightMix (TIB Molbiol)  
27: Commercial assay system / kit \*\*  
28: *In house* PCR assay      29: Other \*\*

#### Group [III] (Detection / identification)

31: Commercial assay system (codes 20-27)  
32: Agarose gel electrophoresis  
33: Hybridization with labelled probe  
34: Nested PCR,      35: Real-Time PCR  
37: DNA sequencing      39: Other \*\*

#### Group [IV] (Positive and/or inhibition control)

41: Commercial assay system (codes 20-27)  
42: Internal control (recombinant plasmid etc.)  
43: External controls (plasmid,  $\beta$ -globin gene, etc.)  
44: Genomic DNA of the corresponding pathogen  
48: Other \*\*      49: No inhibition control applied

Group [V] (Target gene) - - not applicable

#### Group [VI] (PCR results)

61: **Positive**      63: **Questionable**  
62: **Negative**      64: **Inhibition**

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
## BACTERIAL GENOME DETECTION

### PCR-/NAA *C.difficile* toxin genes (RV 545)

Detection of *Clostridium difficile* Toxin-encoding genes by PCR or other procedures for nucleic acid amplification & detection (NAA)

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as culture aliquots or native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *Clostridium difficile* toxin gene PCR/NAA assays established in your routine diagnostic setting. **Note:** The DNA content of the samples may be significantly lower compared to direct extraction or boiling a bacterial colony.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of DNA from *Clostridium difficile* toxin genes in the samples. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to target organisms, thereby allowing inhibition controls targeting human genes like β-globin. Please apply the same methods for DNA extraction and amplification like you are using for the workup of cultures or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of *C. difficile* toxin genes is requested, you are free to report the corresponding toxin genes [codes 70 to 79] in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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### CODE Numbers

PCR-/NAA *C. difficile* toxin genes (RV 545)

Code numbers are subdivided into seven groups. Although providing technical specifications of the applied assays is highly desired for statistical purposes (individual data will be kept confidential), completing groups [I] to [V], and [VII], is voluntary.

Completing of group [VII] is mandatory.

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#### Group [II] (Amplification)

20: XPert *C.difficile* (Cepheid)  
21: BD MAX Cdiff (BD)  
27: Commercial assay system / kit \*\*  
28: *In house* PCR assay      29: Other \*\*

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#### Group [V] (Target gene)

51: Commercial assay system (codes 20-27)  
52: *C.difficile* toxin genes *tcdA/ tcdB/ tcdC*  
59: Other\*\*

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optional: 71: *tcdA*      72: *tcdB*  
73: binary toxin      74: *tcdC* (wildtype)  
75: *tcdC* $\Delta$ 117 (deletion)      79: Other\*\*

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

## BACTERIAL GENOME DETECTION

### PCR-/NAA VRE (Vancomycin restant enterococci) (RV 546)

Detection of vancomycin resistant *Enterococcus* spp. DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)  
Instructions for testing and Code numbers (2 pages)

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as culture aliquots or native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and VRE PCR/NAA assays established in your routine diagnostic setting. **Note:** The DNA content of the samples may be significantly lower compared to direct extraction or boiling a bacterial colony.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of vancomycin resistant *Enterococcus* spp. DNA from the samples. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to target organisms, thereby allowing inhibition controls targeting human genes like  $\beta$ -globin.

Please apply the same methods for DNA extraction and amplification like you are using for the workup of cultures or routine clinical samples. Protocol and workflow characteristics should be specified in the Report Form according to the provided code numbers. Although only reporting of qualitative results for the presence or absence of VRE DNA is requested, you are free to report the corresponding *van*-genes [codes 70 to 79] in the Comments column. The latter information is voluntary, but very helpful for the final evaluation of this proficiency testing program.

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## FUNGAL GENOME DETECTION

### PCR-/NAA *Pneumocystis jirovecii* (RV 560)

Detection of *Pneumocystis jirovecii* DNA by PCR or other procedures for nucleic acid amplification & detection (NAA)

Instructions for testing and  
Code numbers (2 pages)

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#### Specimen preparation

Samples should be stored at 4°C upon receipt. The closed tubes should be briefly centrifuged to pellet the lyophilized material. Immediately before testing, renature the samples by adding **300 µl** of sterile water (PCR-grade) and incubate at room temperature for 20 min. on an orbital shaker and/or occasional vortexing. When the pellets are completely dissolved, the resulting suspensions should be considered as native clinical specimens and 100 µl aliquots processed as such using typical protocols for DNA extraction and *P. jirovecii* - PCR/NAA established in your routine diagnostic setting.

#### Testing and reporting of results

The concept of this proficiency panel is designed for determination of *P. jirovecii* DNA in the sample material. The sophisticated composition of the sample matrix comprises of proteins, bacteria, and human cells next to the target organisms, thereby allowing for the application of optional inhibition controls targeting human genes (like  $\beta$ -globin, etc.).

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#### Group [II] (Amplification)

20: MycAssay Pneumocystis (Myconostica)  
21: LightMix *Pneumocystis jirovecii* (TIB Molbiol)  
22: RIDAGENE *P. jirovecii* (r-Biopharm)  
23: AmpliGnost *P. jirovecii* PCR Kit  
24: *P. jirovecii* Real TM (Sacace)  
27: Commercial assay system / kit \*\*  
28: *In house* PCR assay      29: Other \*\*

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31: Commercial assay (please specify)  
32: Agarose gel electrophoresis  
33: Hybridization with labelled probe  
34: Nested PCR  
35: Real-Time PCR (TaqMan format)  
36: Real-Time PCR (LightCycler)  
37: DNA sequencing      39: Other \*\*

#### Group [IV] (Positive and/or inhibition controls)

41: Commercial assay (please specify)  
42: Internal control (recombinant plasmid etc.)  
43: External controls (plasmid,  $\beta$ -globin gene, etc.)  
44: Genomic DNA of the corresponding pathogen  
48: Other \*\*      49: No inhibition control applied

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